

## ChEC in budding yeast

Gabe Zentner

Modified from Schmid et al, *Mol Cell*, 2004

---

1. Grow 50 mL cells to  $OD_{600} = 0.5-0.7$  in a 300 mL flask
  2. Harvest cells in a 50 mL conical tube at  $1,500 \times g$  for 1 min
  3. Wash cells 3 x 1 mL Buffer A. Transfer cells to a 1.5 mL tube with the first wash and spin as above between washes
  4. Permeabilize cells. Resuspend pellet in 600  $\mu$ L Buffer A + 0.1% digitonin (add 30  $\mu$ L 2% digitonin in DMSO to 570  $\mu$ L Buffer A) and incubate at 30°C for 5 min. Remove 100  $\mu$ L as zero timepoint prior to step 5
  5. Add 1.1  $\mu$ L 1 M  $CaCl_2$  (~2 mM final), mix, and incubate at 30°C
  6. At each desired time point, remove a 100  $\mu$ L aliquot of the digestion to a tube containing 90  $\mu$ L stop solution and 10  $\mu$ L 10% SDS and vortex to mix
  7. Add 80  $\mu$ g proteinase K. Digest protein at 55°C for 30 min
  8. Extract nucleic acids. Add 200  $\mu$ L phenol/chloroform/isoamyl alcohol, mix well, and spin at max. speed for 5 min in a microfuge. Transfer aqueous phases to new tubes, add 2  $\mu$ L (30  $\mu$ g) Glycoblue and 500  $\mu$ L 100% ethanol, mix, and precipitate at -80°C for  $\geq 30$  min
  9. Spin at max speed and 4°C for 10 min
  10. Wash pellets with 1 mL 100% ethanol and aspirate ethanol
  11. Briefly air-dry pellets and resuspend in 30  $\mu$ L Qiagen EB or 0.1X TE buffer, pH 8.0
  12. Add 10  $\mu$ g RNase A and incubate at 37°C for 20 min. Run 5  $\mu$ L RNase-treated DNA on a 1.5% agarose gel to check DNA fragmentation if desired
  13. Enrich small DNA fragments. Add 75  $\mu$ L Ampure XP beads to RNase-treated DNA and pipet up and down 10X to mix. Incubate 5 min at RT and collect beads on magnetic rack. Reserve supernatant (unbound fraction, ~105  $\mu$ L, containing smaller DNA fragments) to a new tube containing 96  $\mu$ L Qiagen EB or 0.1X TE buffer and 4  $\mu$ L 5 M NaCl
  14. Extract DNA from the unbound fraction. Add 200  $\mu$ L phenol/chloroform/isoamyl alcohol, mix well, and spin at max. speed for 5 min in a microfuge. Transfer aqueous phases to new tubes, add 2  $\mu$ L (30  $\mu$ g) Glycoblue and 500  $\mu$ L 100% ethanol, mix, and precipitate at -80°C for  $\geq 30$  min
  15. Spin at max speed and 4°C for 10 min
  16. Wash pellets with 1 mL 100% ethanol and remove ethanol with vacuum
  17. Briefly air-dry pellets and resuspend in 25  $\mu$ L 0.1X TE buffer, pH 8.0 or Qiagen EB (10 mM Tris pH 8)
- 

## Solutions

2% digitonin

Add 20 mg high-purity digitonin to 1 mL DMSO. Vortex ~30 sec to dissolve and store 100  $\mu$ L aliquots at -20°C.

Buffer A (100 mL)

1.5 mL 1 M Tris pH 7.5	-	15 mM Tris
8 mL 1 M KCl	-	80 mM KCl
50 $\mu$ L 0.2 M EGTA	-	0.1 mM EGTA

For 10 mL Buffer A (2 samples), add the following before use:

10  $\mu$ L 10 mg/mL LPC (10  $\mu$ g/mL final) or 1 cOmplete EDTA-free mini tablet  
100  $\mu$ L 100 mM PMSF (1 mM final)  
10  $\mu$ L 200 mM spermine (0.2 mM final)  
5  $\mu$ L 1 M spermidine (0.5 mM final)

2X Stop buffer (100 mL)

8 mL 5 M NaCl	-	400 mM NaCl
4 mL 0.5 M EDTA	-	20 mM EDTA
2 mL 0.2 M EGTA	-	4 mM EGTA